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
*Every dentist's dream come true:
always the right colour at hand!"*

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
Amalgam tattoos and endocavities can be reliably covered with the Blocker. The layer thickness depends on the degree of discolouration and the depth of the cavity."

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
"The addition of a flow variant to this product idea is absolutely to be welcomed and rounds off the application possibilities excellently."

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The 9-year-old patient Y. came to our practice after a trauma in the anterior region. It was determined that it was a simple accident (the boy had fallen while playing with friends).




The subsequent examination revealed a fracture (Ellis Class III with pulp exposure) of tooth 1.1, physiological mobility of the tooth and other teeth, mild swelling of the upper lip, mild pain and poor plaque management. (Fig. 1-2) The patient had




brought the broken fragment, which was immediately placed in physiological saline solution.

As the pulp exposure was minor and not bleeding, it was decided to fix the fragment again and reconstruct the missing part directly.

After local anaesthesia with articaine, the surgical area was prepared. Placement of the rubber dam was quite difficult due to the mixed dentition.

The surface was carefully cleaned with a toothbrush. The adhesive system Tokuyama EE-Bond (7th generation with enamel etching) was then applied to the tooth and the fragment. The fragment was repositioned with a thin interposition of flowable composite. (Figs. 3-5)

After curing (Fig. 6), a bevel was created along the fracture line vestibularly and palatally to increase the adhesion of the fragment (Figs. 8-9), smooth the edges and allow direct reconstruction of the missing part.

After reapplication of the adhesive system (Fig. 10), a freehand layering was carried out. (Fig. 11)

The first aesthetic examination had revealed a clear incisal translucency and a distinct halo at the edge.

NE (Tokuyama Asteria) was applied as the first palatal layer. Afterwards, the build-up was done with Body A2.

According to the manufacturer's instructions, Body was applied up to the middle third, leaving space for a thin layer of enamel NE on the incisal third. (Fig. 12)

To create the enamel translucency, TE was used between the body incisal cusps. The halo was created with Body A2. (Fig. 13)

After finishing (Fig. 14) and polishing (Fig. 15), the patient was discharged home and regularly questioned about possible pulpal symptoms in the following days.

The control image after 7 days (Fig.

16) shows the excellent integration of the restoration. Due to the integration of the composite used and the simple layering technique, we were able to achieve a satisfactory result immediately, so that no further treatment was required.

The body shade did an excellent job of covering the transition line between the tooth and the restoration. With the help of enamel-NE and translucent TE, an absolutely natural incisal area has been created.

The patient was satisfied with the aesthetic result of the restoration. The last picture shows the postoperative check after one week and rehydration of the adjacent teeth.

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
Customers also bought this additional product:
UNIVERSAL BOND II

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

Customers also bought this additional product:
BOND FORCE II

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
Always the right viscosity with perfect polishability make this flow the agent of choice."

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


Shades - 5:

Customers also bought this additional product:
BOND FORCE II PEN

[More at BONDINGS >> BOND FORCE II PEN](#)

Tooth forming with composite

Tooth reshaping with composite after tooth transplantation in the anterior region Prof. Dr. Anne-Katrin Lühns, Hannover

History: After a fall in a swimming pool in 2000 and a total dislocation of the two central incisors 11 and 21 (both teeth were lost in the fall), teeth 34 and 44 were transplanted in region 11 and 21. Seven years later, the patient, who had been satisfied with the restoration up to that point, presented for a consultation regarding the further restoration of the transplanted premolars. Orthodontic treatment with regular radiographic control of region 11/21 was carried out in parallel. The apical region of the transplanted premolars was inconspicuous, there was no evidence of root resorption.

Treatment planning: After evaluating the clinical and radiographic findings, the plan was to minimally invasively reshape the grafted premolars 34 and 44 into central incisors. Other treatment alternatives such as the fabrication of all-ceramic veneers and the performance of a surgical crown lengthening (vestibular) to harmonise the gingival line were rejected by the patient for reasons of invasiveness.

Treatment procedure: First, a wax-up (Fig. 2) was made to inform the patient about different treatment options and to evaluate the expected loss of substance due to preparation measures.

Based on the silicone key fabricated on the wax-up, an intraoral mock-up

to shape the restoration more delicately in a vestibular direction.

The teeth were now cleaned with fluoride-free prophylaxis paste. Then the shade was selected and the material Estelite Sigma Quick (Tokuyama, Altenberge) was used as the composite system.

Fig. 1a, b: Transplanted premolars in region 11/21, composite copings from orthodontic treatment still in situ

Fig. 2: Wax-up of region 11/21 and drawn-in suggested course of the gingival margin after surgical crown lengthening

was fabricated with a provisional plastic material to visualize the targeted treatment results. Due to the thickness of the wax-up and the position of the transplanted premolars in the dental arch, it was agreed with the patient, who was very satisfied with the mock-up,

Fig. 3a, b: Situation after removal of the composite flaps (a), intraoral mock-up for visualisation of the possible treatment result (b).

This system contains an extra-opaque dentin material, OPA2, which was used to build up the incisal dentin core that is not supported by tooth structure.

Fig. 4a-m: Premolar remodeling T1 (analogous procedure was followed for T2): wax-up with silicone key (a), silicone key in situ (b), enamel conditioning (c), etching pattern vestibular (d) and palatal (e), flowable in the silicone key (f), built-up posterior wall (g), layering of the dentine core (h, i), enamel mass T1 (j) and T2 (k), restorations after final polishing of vestibular (l) and palatal (m).

Fig. 5a, b: Clinical situation 4 weeks postoperatively, view from vestibular (a) and palatal (b)

Fig. 6a, b: Reshaped premolars after 4.5 months postoperatively, view from vestibular (a) and palatal (b)

Fig. 7a, b: Composite restorations after 14 months postoperatively, vestibular view (a) and palatal view (b)

The following shades were used: OPA2, OA2, OA3 (dentin shades), A2, A3 (enamel masses), CE (Clear Enamel, effect shade for translucent effects), OPA2 (Flowable).

The enamel surface of the grafted premolars was minimally invasively roughened with a finishing diamond and then conditioned for 30 seconds using 37% orthophosphoric acid (Conditioner 36, Dentsply, Constance) and adhesively pretreated (Heliobond, Ivoclar Vivadent, Ellwangen).

With the help of a silicone key, the posterior wall was now built up from the flowable of shade OPA2 in a thin

layer thickness. The different dentin materials were then modelled in several layers: OPA2 as opaque core, OA3 cervical and OA2 incisal.

An approx. 1 mm wide area was left free in the area of the incisal edge, which was filled with a thin layer of effect material (CE). A very thin layer of OA2 was added directly to the incisal edge to achieve the "halo effect". The dentine core was then coated with the enamel materials (A2 incisal, A3 cervical).

After final polishing, the teeth were fluoridated (Biflourid, VOCO, Cuxhaven).

The illustrations above show the restorations at follow-up appointments after 4 weeks, 45 months and 14 months.

Follow-up sessions were arranged with the patient at 6-month intervals after completion of treatment.

Epicrisis and discussion: The patient case presented represents the reshaping of two premolars after autogenous tooth transplantation following anterior trauma, in which both central incisors were totally luxated and were not available for replantation because they could no longer be found.

With regard to the reshaping of transplanted premolars into anterior teeth, there are individual case reports in which both composite restorations and all-ceramic veneers were fabricated as a therapeutic measure for tooth reshaping.

In this case, minimally invasive composite restorations were fabricated at the request of the patient, who was informed about a wide variety of procedures. The posterior wall was built up from a very thin layer of Flowable using a silicone key after a previous wax-up. Since the occlusal contacts are in the area of the tooth structure, this procedure proved to be practicable.

The patient is satisfied with the treatment result achieved. During the observation period of 14 months, there were no indications of restoration failure.

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
The image displays a collection of Tokuyama dental bonding products. In the foreground, a white box for 'Tokuyama Universal Bond II Self-cured Dental Universal Adhesive KIT' is visible. Next to it are two bottles: a pink one labeled 'BOND A' and a blue one labeled 'BOND B'. To the right, there is a syringe and a pen applicator, both labeled 'BOND FORCE II'. The background features a red box and a blue box, also labeled 'BOND FORCE II'. The entire scene is set against a white background with blue molecular structure graphics.

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
Dentist Gerhard Dalheimer has been in private practice in Schwalbach/Elm for 28 years. Over time, he has specialised in the areas of Cerec and CMD + prosthetics and treats his patients together with his colleague Dr Hasan Ackay.

Efficiency is a very important factor in everyday dental practice and, against the background of constantly rising

costs, materials such as Universal Bond II and ESTECM II PLUS from Tokuyama Dental are good alternatives to previous materials because they enable a universal and correspondingly time-saving method of working. The following four patient cases show how dentist Gerhard Dalheimer used these materials in clinical trials to the complete satisfaction of the practitioner and the patients.

Case 1: Bonding a veneer to tooth 13

The 70-year-old patient presented because the veneer on 13 had come loose (Fig. 1a). After inspection of the tooth and the veneer, the reattachment could be started. The veneer and the tooth were cleaned of the old cement and then it was tried on again.

The veneer and the tooth were then cleaned with alcohol and dried, the Universal Bond II (A/B) was mixed and applied to the veneer and the tooth and blown (Figs. 1b and 1c). The Universal Paste was applied to the veneer and placed on the tooth (Fig. 1d).

The excess was removed after 3 sec. of light curing and then cured again from the vestibular and palatal sides for 1 min. each (Fig. 1e). The result is very satisfactory after a very short treatment (Fig. 1f).

Case 2: Cementation of a post in preparation for a bridge abutment on tooth 14

The 56-year-old patient has been under our care for a long time. She is healthy and her medical history is unremarkable. Tooth 14 has been endodontically restored for a long time and had no prosthetic restoration. The patient presented for a new restoration after a fracture of the old build-up filling (Fig. 2a). We took a radiograph (Fig. 2b) and discussed the treatment options of a core build-up and a bridge restoration. The two root canals were prepared with the Tokudrill up to the first mark (Fig. 2c) and then the length was checked with the Verifier.

Subsequently, the quartz fibre posts were inserted on a trial basis (Fig. 2d and 2e). The Universal Bond II (A/B) was mixed and applied to the pins as well as into the canals and blown (Fig. 2f). The Universal Paste was applied to the canals (Fig. 2g) and the posts were inserted into the canals (Fig. 2h).

The excess was removed and then light-cured. Finally, the die was built up with a dual-curing composite and the tooth and prepared (Fig. 2i to 2l). The bridge restoration will be done shortly.


Case 3: Cementation crown

This patient is 45 years old. Tooth 45 had already been treated endodontically some time ago. In the course of this clinical trial, it has now been restored with a crown in one-single-visit chairside. The crown is made of a zirconium oxide-reinforced lithium silicate ceramic (centra dou). After the tooth was prepared (Fig. 3a), the crown was fitted. First, the crown and the stump were cleaned with alcohol and dried, the Universal Bond II (A/B) was mixed and applied to the crown as well as the stump and blown (figs. 3b and 3c). The Universal Paste was applied to the crown and used on the die (Fig. 3d).

The excess (Fig. 3e) was removed after a 3 sec. light cure. The crown was then cured again from all sides for 30 seconds each. The result after such a short treatment period is very satisfactory for the patient (Fig. 3f).

4a

4b

4c

4d

Case 4: Repair of veneer 13, 12

The patient is 65 years old and presented to our practice because the veneer on the crowns region 12,13 had chipped off (Fig. 4a). First, Universal Bond II (A/B) was mixed (Figs. 4b and 4c), applied to the crowns and air-blown.

Then Estelite Color medium chroma opaque was applied to the non-precious metal framework, followed by restoration of the veneer with Estelite Asteria A2B and finally polishing.

Even such indications, which are very unattractive for our patients from an aesthetic point of view, can be easily solved in a very short time and to the complete satisfaction of the practitioner and patient (Fig. 4d) with the universal adhesive system Universal Bond II from Tokuyama and the other, optimally matched product range from Tokuyama Dental.

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
Tokuyama Estesil is an excellent impression material for daily use"

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FI | UNIDENT AB
& SERVICE LTD
TOKUYAMA Dental Deutschland GmbH